



# Warm and Cold Compression Therapy Machine VU-COT01-V2

Intelligent temperature control **Rapid** cooling/heating



Continuous use for **240 hours**  
Safety and stability



With precooling function  
One click Start



**Iceless**,  
Semiconductor  
Cooling/heating  
technology



2 separate tanks, rapid  
toggle cold/ warm therapy



LCD display full  
touch screen control



## Cold/ Heat Compression Therapy Machine VU-COT01-V2

VU-COT01-V2 This new cold and hot therapy system integrates iceless cold, heat and compression therapies in one easy-to-use unit. The cold compression machine adopts a fully automatic semiconductor refrigeration system, with optional wrapping accessories for different body parts. For patients with long-term exercise, labor, various trauma, bone injury, fever, pain and swelling etc. it can continuously improve the blood of the human body circulation and vasodilation. It has the advantages of safety, simplicity, and high efficiency.

### Specification

Working Voltage: DC12V  
Power Consumption: 400W  
Adjustable Time: 5~90minutes  
Pressure range: 30~120mmHg  
Water temperature:  
3°C~15°C (37°F~59°F) ,  
38°C~40°C (100°F~104°F)  
Machine size: 50\*28\*28.5cm  
Machine net weight: 10KG  
Switchable temperature unit °C and °F



### Component

- Host (including air pump)
- Connecting pipe
- Power adapter
- Sleeves
- trolley case (60x33x42cm)



### Optional 5 types of therapies

 Cold Compress    Heat Compress    Compression    Cold Compression    Heat Compression

### Features & Details



LCD display full  
touch screen control




Water compress 3°C~15°C  
can enjoy iceless cold



With power adapter, can be connected  
to power supply 100V~240V

## Optional accessories for different body parts




## Product's intended use

The equipment is used for postoperative and post-injury rehabilitation treatment, to alleviate inflammatory reaction, control swelling and relieve pain by applying pressure and providing cold and hot compress treatment.

## Market scenario

Market scenario: It can be used in clinical departments such as orthopedics, sports medicine, and rehabilitation, as well as team doctors of various sports teams.



XIAMEN WEIYOU INTELLIGENT TECHNOLOGY CO.,LTD



TEL:(86)-592-6251545

FAX:(86)-592-6252548

ADD:Unit 301, No. 6 Xianghong Road,Torch Hi-Tech Zone Industrial Park,Xiang'an District,Xiamen P.R.China.

WEB:www.weiuit.com  
MAI:sales@weiuit.com



# Warm and Cold Compression Therapy Machine VU-COT01-V2